

# Using NovoCyte® Benchtop Flow Cytometer for Multiplexing Antibiotic Screening in Food Commodities with BEADYPLEX™ Assay

Nan Li<sup>1</sup>, Xiaobo Wang<sup>1</sup>, Celia Suárez Pantaleón<sup>2</sup>, Alessandro Gallo<sup>2</sup>, Benoit Granier<sup>2</sup>  
<sup>1</sup>R&D, ACEA Biosciences, Inc., San Diego, CA, United States (nli@aceabio.com)  
<sup>2</sup>Unisensor S.A., Liège (Ougrée) - BELGIUM (celia.suarez@unisensor.be)

## Introduction

The misuse of antibiotics in animal farming may lead to the presence of residues in edible products, thus compromising consumer health (toxicological and allergic reactions, as well as bacterial antibiotic resistance development) and industrial food processes. Rapid and efficient screening analytical methods for the early detection of antibiotic residues are therefore essential to guarantee food safety and to reduce the negative impact on industrial fermentation processes.

The current trend in food analysis entails the implementation of multi-residue technologies which considerably improve analysis time and overall cost. BEADYPLEX™ is a competitive suspension bead-based Flow Cytometric ImmunoAssay (FCIA) for the simultaneous screening of more than 80 residues from the 10 most relevant antibiotic families, namely tetracyclines, sulfonamides, β-lactams, aminoglycosides, macrolides, fluoroquinolones, lincosamides, phenicols, polymyxins and pleuromutlins, providing the family identification in one single test. The method has been validated in muscle (porcine, bovine and poultry) and fish (salmon and coley) following the 2010 CRL Guidelines supplementing European Decision 2002/657/EC. Preliminary results show the potential applicability of BEADYPLEX™ to other matrices, such as milk or eggs.

BEADYPLEX™ has been validated using NovoCyte® 2000 flow cytometer (with 488 nm and 640 nm lasers and 4 fluorescence channels) and the BeadyExpress™ software developed by ACEA Biosciences. NovoCyte® is a high performance bench flow cytometer featured with wide dynamic range (7 logs), high sensitivity, high speed, and ease of use. BeadyExpress™ software, developed for BEADYPLEX™ FCIA multiplexing assay, allows the straightforward and automatic data acquisition and analysis. With a few clicks of button, the software automatically sets the gates for each population of the BEADYPLEX™ beads, acquires and analyses the external fluorescence signal, and confirms the presence/absence of antibiotics, using the proprietary developed algorithm. The software generates a report indicating the validity of the assay and identifying each analyzed sample as Positive or Negative for the 10 antibiotic families.

## Assay Principle and Methodology



A twelve-plex competitive Flow Cytometry Immunoassay has been developed for the detection of 10 antibiotic families in food commodities (more than 80 residues included).

Antibiotic-mimicking molecules are individually coated on fluorescent beads encoded by their size and internal fluorescence. First, a mixture of primary binders directed to the different antibiotic families are incubated with the sample extracts and the antibiotic-conjugated beads. Then, secondary binders labeled with a fluorescent reporter are added for incubation. A washing step is required after every incubation step.

Upon reading with the NovoCyte®, the beads are classified according to their size and internal fluorescence (excited by 640 nm laser and detected in APC channel). The external assay signal for every bead population is detected from the secondary binders in FITC channel by 488 nm laser excitation.

According to the competitive principle of the assay, the absence of antibiotics in the sample is translated into the maximum fluorescent assay signals, whereas a decrease in fluorescence is observed upon the presence of free antibiotics for the associated bead.

### Flow Cytometry Immunoassay



### Sample Extraction

1 g sample + 1 mL extraction buffer (10 min) → centrifugation + filtration

### 1st Incubation Step

Sample extract + antibiotic-coated beads mix + primary antibodies mix 30 min, RT, darkness

### Washing

### 2nd Incubation Step

Fluorescent secondary reporters 15 min, RT, darkness

### Washing



### Flow Cytometer Reading

NovoCyte® + BeadyExpress™

## Results

Most of the antibiotics included on the scope of BEADYPLEX™ are detected at concentrations equal to or below the European regulatory limits. LOD (Limit of Detection, Detection Capability).

Family	Antibiotic	LOD (µg/kg)	Family	Antibiotic	LOD (µg/kg)	Family	Antibiotic	LOD (µg/kg)												
Aminoglycosides	Streptomycin	250	Tetracyclines	Sulfadiazine	50	β-Lactams	Chlortetracycline	25												
	Dihydrostreptomycin	125		Sulfamethoxazole	10		Doxycycline	25												
	Gentamicin	50		Sulfamethoxazole	10		Oxytetracycline	50												
	Neomycin B	100		Sulfadimethoxine	50		Tetracycline	50												
	Kanamycin A	100		Sulfadoxine	10		Demeclocycline	50												
	Paromomycin	250		Sulfachloropyridazine	1		Methacyclocline	25												
	Sulfonamides	Cefquinome		50	Quinolones		Marbofloxacin	20	Macrolides	Flumequine	800									
		Ceftiofur		50			Ciprofloxacin	0.5		Erofloxacine	5									
		Desfuroylceftiofur		2000			Sulfafloxacin	0.5		Danofloxacin	250									
		Cefoperazone		3			Sulfamethoxyypyridazine	0.25		Oxolinic acid	800									
Cefalexin		3500	Sulfamerazine	20		Difloxacin	50													
Cefalonium		5	Sulfamethazole	1		Norfloxacin	2.5													
Cefepim		75	Sulfamonomethoxine	10		Saralofloxacin	125													
Desacetylcefepim		100	Sulfquinoloxaline	50		Ofloxacin	10													
Cefazolin		100	Sulfafurazole	5		Enoxacin	50													
Cefazolin		50	Phenicols	Pleuromutlins		Polymyxin	Colistin	Colistin		250										
Penicillin V	12.5	Phenicols			Pleuromutlins				Polymyxin		Colistin	Colistin	250							
β-Lactams	Penicillin G													Ampicillin	Amoxicillin	Oxacillin	Cloxacillin	Dicloxacillin	Nafcillin	Piperacillin

Example of automatically generated report by the software BeadyExpress™ upon plate reading with NovoCyte®.



### Blank and Positive Control Samples

Validation of the test

### Unknown Samples

Identification of Positive **POS** and Negative **NEG** results for each antibiotic family

## Conclusions and Perspectives

BEADYPLEX™ is a user-friendly assay that provides sensitive and specific broad spectrum detection of antibiotics in multiple matrices in one single test with high-throughput capability.

NovoCyte® and BeadyExpress™ software highly simplify the data acquisition and the interpretation of results, by providing a direct identification of the antibiotic family from the sample, with significant reduction of analysis time and cost.

The work presented here illustrates the potential of using Flow Cytometry for food analysis. New FCIA kits are under development in order to expand the application of this versatile and flexible technology to other types of contaminants.

**Acknowledgement** • This work has received funding from Eureka and the Belgian Walloon Region. Beadyplex™ kits have been developed with the CER Groupe (Marloie, Belgium).



ACEA Biosciences, Inc.  
 6779 Mesa Ridge Road #100  
 San Diego, CA 92121 USA  
 Tel: +1 866-308-2232  
 Fax: +1 858-724-0927  
 info@aceabio.com

Allée de la Cense-Rouge, 98  
 4102 Ougrée (Liège) - Belgium

+32 4 252 66 02  
 +32 4 252 90 55

info@unisensor.be  
 www.unisensor.be